

Anuform®

Intra-Anal and small Intra-Vaginal Probe


The revolutionary contoured shape of the Anuform® provides a sound electrode/tissue interface, which allows for comfortable muscle stimulation and accurate biofeedback. As the Anuform always sits in the same position, comparable EMR readings can be taken over different sessions.



How does the Anuform® work?

When the Anuform is inserted, the anal sphincter closes around the narrow neck of the probe. This secures the correct positioning without the need to hold the probe in place.

- Revolutionary shape for repeatable positioning and ambulatory use
- Can be used with most muscle stimulators and biofeedback equipment (refer to manufacturers instructions before use)
- Single patient use
- Single orifice use
- Medical grade surgical stainless steel



Product Code 091253640